GB Instructions for use/Technical description
Caiman 5 Articulating, Sealing & Cutting Instruments

USA Instructions for use/Technical description
Caiman 5 Articulating, Sealing & Cutting Instruments
Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
Aesculap® Caiman 5 Articulating, Sealing & Cutting Instruments

Legend

Instrument

- 1 HF activation button
- 2 Rotation wheel
- 3 Jaw articulation
- 4 Moveable jaw part with marking
- 5 Articulation activation lever
- 6 Blade activation lever
- 7 Jaw part activation lever
- 8 HF cable with connector

Generator

- 9 SEE LEKTRAFUSE HF GENERATOR IFU
- 10 Foot switch connection socket
- 11 Regrasp signal lamp
- 12 Display
- 13 Error signal lamp
- 14 Instrument connection socket
- 15 "Power on" indicator light
- 16 Foot control switch

Symbols on product and packages

	Sterilization with irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product
	Follow the instructions for use
	Date of manufacture

Scope

Aesculap Caiman 5 Articulating, Sealing & Cutting Instruments are sterile, disposable instruments. Caiman 5 instruments with the Aesculap Lektrafuse HF generator form an HF system for applying HF energy. Caiman 5 instruments can only be used in conjunction with the Lektrafuse HF generator.

The user can seal and cut tissue to a length of up to 26.5 mm with Caiman 5 instruments. The energy input into the tissue is controlled by the Lektrafuse generator software.

- ▶ For item-specific instructions for use and information regarding material compatibility, see Aesculap Extranet at www.extranet.bb Braun.com

Intended use

Indications

Caiman 5 instruments are used for grasping, dissecting, sealing and cutting tissue during open and laparoscopic surgical procedures.

Caiman 5 instruments can be used on vessels with a diameter of up to and including 7 mm.

Contraindications

- There is no evidence that the Caiman 5 instruments are suitable for tube sterilization, coagulation of the fallopian tubes, etc. Therefore the product is strictly not to be used in sterilization procedures.
- Use on vessels with a diameter of greater than 7 mm. Particular caution must be exercised in case of pathological changes of vessels with a diameter of less than 7 mm, e.g. atherosclerosis.

As the safe and effective use of electrosurgical instruments is in principle greatly dependent on influencing factors that only the user can control, the above information can only be used as general guidance only. The clinically successful use of such products is dependent on the knowledge and experience of the surgeon, who is responsible for deciding which structures can be reasonably treated, taking into account the safety and warning instructions given in these instructions for use.

Safe handling and preparation

Risk to the patient and/or user when using the Caiman 5 instruments with HF sources other than the Lektrafuse HF generator!
Non-observance can result in injury or death and restriction in the function of the Caiman 5 instruments.

- ▶ Use Caiman 5 instruments only with the Lektrafuse HF generator.

Risk of death or injury to the patient through non-observance of the following directions!

- ▶ Use Caiman 5 instruments with caution on patients with cardiac pacemakers or defibrillators.

Risk of injury due to evaporation of tissue fluid!

During the sealing process, tissue fluid can be converted into steam by the HF energy. This can result in collateral damage to the directly neighboring tissue.

- ▶ Observe the sealing process carefully.

Risk of injury caused by incorrect operation of the product!

- ▶ Do not activate the HF current before the Caiman 5 instrument has been completely closed and the lock is activated.
- ▶ Do not open Caiman 5 instruments during the HF process as this will automatically interrupt the sealing process.

Risk of injury from ignition or explosion of flammable gases!

Sparks may occur when using the Lektrafuse HF generator as directed.

- ▶ Follow the safety notes in the instructions for use of the Lektrafuse HF generator.

Risk of insufficient sealing due to contaminated electrode surfaces!

- ▶ Keep the product's contact surfaces clean during surgery. Remove encrusted tissue residues or body fluids with a moistened swab.
- ▶ When cleaning the electrodes:
 - Do not activate the HF current.
 - Do not immerse the device jaws in liquid, including saline solution.
 - Do not clean the device jaws with abrasives.

Risk of injury to patient or damage to the instrument when removing the Caiman 5 instrument with distally articulated jaw part!

- ▶ Check the compatibility of the trocar and the Caiman 5 instrument.
- ▶ Before inserting or removing the Caiman 5 instrument, place the articulation mechanism in the neutral position.
- ▶ Insert the Caiman 5 instrument into the trocar and remove it only when locked, and carefully.

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite knowledge and experience.

- ▶ Read, follow, and keep the instructions for use.
- ▶ Only use product after consulting the instructions for use of the Lektrafuse HF generator GN200.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.

The product is gamma-sterilized and supplied in sterile packaging.

The product must not be reused.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.

The accessory voltage rating of the product is 240 Vp.

To avoid HF burns:

- ▶ Always keep the working end of the product in the user's field of vision whenever the HF power is activated.
- ▶ Prior to activating the Lektrafuse HF generator, check that the working end of the product is not touching any electrically conductive accessories.
- ▶ Never place the product on or next to the patient.

Safe operation

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.
- ▶ Before activating product, ensure that no electrically conductive objects are in contact.
- ▶ Do not place the product on or next to the patient.

Risk of damage to the tissue or the instrument due to incorrect use of the instrument with the lock engaged!

- ▶ Avoid the following actions when the lock is activated:
 - Rotation of the shaft
 - Articulation of the jaw parts

Risk of insufficient sealing or transection of the tissue!

- ▶ Do not bundle or fold the tissue which is to be sealed and cut, or clamp more than one layer between the jaw parts.
- ▶ Grip the tissue that is to be sealed centrally between the electrodes.
- ▶ Do not start the process if the instrument jaws are in conductive liquids (e.g. blood or saline solution).
- ▶ Do not start the process if there are conductive objects (e.g. clamps, vessel clamps, clips, etc.) between the instrument jaws.

The software of the Lektrafuse HF generator recognizes the Caiman 5 instruments and applies the corresponding device settings. The duration of the sealing process is controlled by the generator software.

The sealing process is started by tapping on an HF button once. The sealing process is interrupted by again tapping an HF activation button.

Note

In case of unsuccessful or incomplete sealing, as well as insufficient functionality of the instruments, a visual and acoustic alarm will be generated by the generator.

Note

In the event of visibly insufficient sealing of the tissue, the sealing quality can be improved by again activating the HF process.

Tissue manipulation, sealing of vessels and tissue bundles and tissue separation

- ▶ Remove the Caiman 5 instrument from the sterile packaging.
- ▶ Insert the plug 8 of the Caiman 5 product into the connecting socket 14 of the Lektrafuse HF generator.
- ▶ In order to bring the jaw parts into the desired position, turn the rotation wheel 2. Whilst doing this, ensure that the black line of the rotation wheel 2 indicates the central position.
- ▶ Bring the tip of the instrument into the desired position for surgery by rotating the shaft and applying an articulation movement to the articulation activation lever 5.
- ▶ Grip the tissue that is to be sealed and/or cut between the jaw parts.
- ▶ Lock the moveable jaw part 4 using the activation lever 7.
- ▶ In order to seal the clamped tissue, actuate the HF activation button 1 on the instrument or foot switch 16.

Note

For additional options see the instructions for use of the Lektrafuse HF generator GN200.

- ▶ In order to cut clamped and sealed tissue, use the blade activation lever 6.

Risk to patient and damage to the instrument due to clamping of metallic objects!
▶ Do not seal or cut over clips or clamps.

- ▶ If an error report appears:
 - Do not cut clamped tissue.
 - Loosen the lock of the handle and open the jaw part.
 - When sealing quality is sufficient, the sealed tissue can be cut using surgical scissors or the instrument.
 - When sealing quality is insufficient, the sealing process must be restarted. If the additional HF cycle leads to a normal seal, the tissue can be cut, and the instrument can be opened and removed.
- ▶ Do not leave the jaw actuation lever in the locked position when the instrument is not in use.

Single-use products

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!
▶ Do not reprocess the product!

Storage

- ▶ Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Troubleshooting

Errors, faults and warnings are displayed on the Lektrafuse HF generator as follows:

- An individual sound or a clear sound sequence and interruption of the sound sequence of HF operation
- The yellow regrasp error light or the red error light on the front side of the Lektrafuse HF generator lights up
- Notification in the display on the front side of the Lektrafuse HF generator
- ▶ For details on errors, faults, warnings and remedying such, see the instructions for use of the Lektrafuse HF generator.

Technical Service

Risk of injury and/or malfunction!
▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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