
### Aesculap<sup>®</sup> Single Use Surgical Instruments

Procedure Sets Instruments for Ward Use and Emergency Care



Aesculap Surgical Technologies

A new Standard in Surgical Instruments







#### Pack per case – The special set

B. Braun Aesculap procedure sets make a difference, independently if used for small procedures or for treatment of large wounds. Our "Pack per case" range of procedure sets include all the single use products which hospitals, clinics and practitioners rely upon every day. However the surgical instruments in the sets make the difference. They are high performance precision instruments, although designed for single use. No wonder, Aesculap instruments represent quality, economic efficiency and reliability since generations.

#### SUSI<sup>®</sup> – The difference

SUSI® (Single Use Surgical Instruments) Instrumentation is a special range of single use surgical instruments. Utilising proven established design patterns and the use of light high performance polymers for the production provides you with high precision instruments with excellent mechanical properties. The material used allows simple and environmentally friendly incineration utilising existing disposal systems in hospitals.

SUSI® instruments are also available individually sterile packed in a dispenser box.



#### The choice is yours – Establish a new standard

"Pack per case" procedure sets have proven their suitability on hospital wards, in intensive care units, A&E Departments, Outpatient Departments, Nursing homes, GP surgeries & practices and can also be used by ambulance services. For all these areas of application, SUSI® procesure sets not only offer medical advantages but also plenty of savings potential. Should you have specific needs and requirements (i.e custom made sets or economic analysis) please do not hesitate to contact us.

Single Use Surgical Instruments



Contents	Pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Gauze swabs, plumsize	5
SUSI <sup>®</sup> forceps Adson, sharp, 120 mm	1
SUSI <sup>®</sup> needle holder Mayo-Hegar, 150 mm	1
SUSI <sup>®</sup> scissors, sharp/blunt, 145 mm	1
SUSI® haemostatic forceps Pean, straight, 145 mm	1

Order unit: PAK = Pack of 20 sets

Contents	Pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Gauze compress, 7.5 x 7.5 cm, 8-fold (ply)	5
Gauze swabs, plumsize	5
SUSI® forceps Adson, sharp, 120 mm	1
SUSI® needle holder Mayo-Hegar, 150 mm	1
SUSI <sup>®</sup> scissors, sharp/blunt, 145 mm	1

Order unit: PAK = Pack of 20 sets

### 5061020 Surgical Dressing Set 2



# 5061030 Surgical Dressing Set 3

Contents	Pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Drape, 75 x 90 cm	1
Hole towel, 75 x 75 cm, adhesive area	1
Gauze compress, 7.5 x 7.5 cm, 12-fold (ply)	10
Gauze swabs, plumsize	5
SUSI® forceps Adson, sharp, 120 mm	1
SUSI® needle holder Mayo-Hegar, 150 mm	1
SUSI® scissors, sharp/blunt, 145 mm	1
SUSI® haemostatic forceps Pean, straight, 145 mm	1

Order unit: PAK = Pack of 20 sets

Contents	Pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Gauze compress, 7.5 x 7.5 cm, 8-fold (ply)	5
Gauze swabs, plumsize	5
SUSI® tissue forceps, sharp, 145 mm	1
SUSI <sup>®</sup> scissors, sharp/blunt, 145 mm	1

Order unit: PAK = Pack of 20 sets

Single Use Surgical Instruments



Contents	Pieces
Defix kidney tray, white	1
Mini compress, 25 x 25 mm	17
SUSI <sup>®</sup> forceps, sharp, ADSON, 120 mm	1
SUSI <sup>®</sup> scissors, sharp/sharp, 110 mm	1

Order unit: PAK = Pack of 20 sets

Contents	Pieces
Defix kidney tray, white	1
Cutfix Stitch Cutter	1
Gauze compress, 5 x 5 cm, 8-fold (ply)	2
SUSI® dissecting forceps Adson, fine, 120 mm	1

Order unit: PAK = Pack of 20 sets

### 5062000 Thorax Drainage Set, small



Contents	Pieces
Plastic tray, blue	1
Vasco surgical gloves, size 7.5	1
Drape, 45 x 75 cm	1
Hole towel, 75 x 110 cm, blue	1
Gauze compress with slit, 10 x 10 cm, 12-fold (ply)	2
Gauze compress, 10 x 10 cm, 12-fold (ply)	5
Cutfix disposable scalpel size 11	1
SUSI® Metzenbaum dissecting scissors, curved, 180 mm	1
SUSI <sup>®</sup> needle holder Mayo-Hegar, 150 mm	1
SUSI® tubing clamp, 135 mm	1
Suture Premilene® 2/0, 75 cm, DS39 (M)	2
Order unit: PAK = Pack of 5 sets	

Single Use Surgical Instruments

#### 5062010 Thorax Drainage Set, large



#### Plastic tray, 31 x 28.5 x 5 cm 1 Drape, 75 x 90 cm 1 Surgical gloves, size 7.5, powder free 1 Gauze compress, 10 x 10 cm, 12-fold (ply) 8 Slit compress, 10 x 10 cm, 12-fold (ply) 2 Bowl, round, 250 ml 1 SUSI® dressing forceps, curved, 265 mm 1 Drape, 75 x 100 cm, self adhesive 2 Injection syringe, Luer, 10 ml 1 1 Sterican<sup>®</sup> cannula size 1, G 20 x 1 <sup>1</sup>/<sub>2</sub>, Ø 0.90 mm, 40 mm, yellow Sterican<sup>®</sup> cannula, LL, G 18 x 1 <sup>1</sup>/<sub>2</sub>, 1 Ø 1.20 mm, 40 mm, pink Sterican<sup>®</sup> cannula, LL, G 20 x 2 <sup>3</sup>/<sub>4</sub>, 1 Ø 0.90 mm, 70 mm, yellow

Contents	Pieces
Sterican® cannula size 12, G 22 x 1 ¼, Ø 0.70 mm, 30 mm, black	1
Cutfix disposable scalpel size 11	1
SUSI <sup>®</sup> tissue forceps, sharp, 145 mm	1
SUSI® Metzenbaum dissecting scissors, 180 mm	1
SUSI® needle holder Mayo-Hegar, 150 mm	1
SUSI® tubing clamp, 135 mm	1
Suture Premicron green 0 (3.5), 75 cm, DS 30 (M)	1
Suture Dafilon blue 3/0 (2), 75 cm, 1 x DS 19	1
Askina® fix fixation plaster, 15 x 20 cm	2
Order unit: PAK = Pack of 5 sets	

Developed in cooperation with: Dr. med. Thomas Kiefer, Klinikum Offenburg, Ebertplatz 12, 77654 Offenburg, Germany

### 5062020 Thorax Drainage Set



Contents	Pieces
Plastic Tray, light blue, 250 x 137 x 57 mm	1
Hole towel, 75 x 110 cm, blue	1
Askina® fix fixation plaster, 15 x 20 cm	1
Slit compress, 10 x 10 cm, 12-fold (ply)	2
Askina <sup>®</sup> gauze compress, 12-fold (ply), 10 x 10 cm	5
SUSI® sponge forceps, curved, 200 mm	1
SUSI® tissue forceps, 1 x 2 teeth, 145 mm	1
SUSI <sup>®</sup> dissecting scissors Metzenbaum, curved, 180 mm	1
SUSI® needleholder Mayo-Hegar, straight, 150 mm	1
SUSI® tubing clamp, 135 mm	1
Disposable scalpel size 11	1
Single use razor	2
Injection syringe, Luer, 10 ml	1

Contents	Pieces
Sterican® cannula size 1, G 20 x 1 1/2, Ø 0.90 mm, 40 mm yellow	1
Sterican® cannula size 12, G 22 x 1 1/4, Ø 0.70 mm, 30 mm black	1
Sterican° cannula, deep intra-muscular, LL G 20 x 2 3/4, Ø 0.90 mm, 70 mm, yellow	1
Drape 78 x 90 cm	1
Suture Premilene® 2/0 (3) 75 cm, DS39 (M)	2
Suture Dagrofil <sup>®</sup> green 0 (3.5) 75 cm, 1 x DS30	1

Order unit: PAK = Pack of 5 sets

Single Use Surgical Instruments



Contents	Pieces
Struma-Drape, T-Form, 198/259 x 30 cm	1
Tray 28 x 38 x 3 cm	1
Sterican <sup>®</sup> cannula size 2, G 21 x 1 1/2, Ø 0.80 mm, 40 mm, green	1
Sterican <sup>®</sup> cannula size 12, G 22 x 1 1/4, Ø 0.70 mm, 30 mm, black	1
Cutfix disposable scalpel size 10	1
Disposable electrode handle	1
Suture Safil® violet 3/0 (2) 70 cm, HR26	2
Suture Monosyn° undyed 4/0 (1.5) 45 cm, 1 x DSM19	1
Injection syringe, Luer, 10 ml	2
Askina <sup>®</sup> Soft steril, 9 x 10 cm	1
SUSI® dissecting forceps Adson, fine, 120 mm	1
SUSI® tissue forceps, 1 x 2 teeth, 145 mm	2
SUSI® dissecting scissors Metzenbaum, curved, 180 mm	1

Contents	Pieces
SUSI® needleholder, Mayor-Hegar, 150 mm	1
SUSI® self-retaining retractor Weitlaner, 3 x 4 teeth, 130 mm	1
SUSI <sup>®</sup> retractor Roux, 170 mm	1
Gauze compress, 16-fold (ply), 10 x 10 cm, X-ray marker	20
Operation gown, Size L	3
Tablebag, 75 x 141 cm	1
Cover drape, 152 x 193 cm	1

Order unit: PAK = Pack of 5 sets

### 5063010 SUSI<sup>®</sup> Cardiac Pacemaker Set



Contents	Pieces
Cover drape 236 x 300 cm, with slit 9 x 108 cm	1
Anesthesia drape 152 x 250 cm	1
Cover drape, self-adhesive, 75 x 90 cm	2
C-Bend cover 70 + 43 x 200 cm	1
Tray, 33 x 30 x 6 cm, white	1
Bowl round, 250 ml, light blue	1
Gauze swabs, round, 40 x 40 cm, X-ray marker	5
Gauze compress, 12-fold (ply), 10 x 10 cm, X-ray marker	20
Cutfix disposable scalpel size 22	1
Omnifix injection syringe, Luer, Lock, 20 ml	1
Sterican <sup>®</sup> cannula size 1, G 20 x 1 1/2, Ø 0.90 mm, 40 mm, yellow	1
Sucker bag 33 x 41 cm, 2 chambers, self-adhesive	1
Disposable electrode handle	1
SUSI® Mosquito forceps, curved, 125 mm	2
SUSI® dissecting scissors Metzenbaum, curved, 180 mm	1
SUSI® dissecting scissors Metzenbaum, curved, 180 mm	1

Contents	Pieces
SUSI® haemostatic forceps Birkett, curved, 185 mm	2
SUSI® needleholder Crile Wood, straight, 200 mm	1
SUSI® tissue forceps Waughs, 200 mm	1
SUSI® sponge forceps, curved, 265 mm	1
SUSI® dissecting forceps Cushing, 200 mm	1
SUSI <sup>®</sup> self-retaining retractor Weitlaner, 3 x 4 teeth, 130 mm	1
SUSI® retractor Roux, 170 mm	1
Bowl, round, 120 ml, transparent	2
Askina <sup>®</sup> Soft 9 x 10 cm	1
Suture Monosyn® violet 3/0 (2), 70 cm, HR26	1
Suture Monosyn <sup>®</sup> undyed 4/0 (1.5), 70 cm, DSMP24	1
Suture Silkam <sup>®</sup> black 0 (3.5), 75 cm, 1 x HR26	1
Operation gown, Size L	2
Table bag, 75 x 141 cm	1
Table cover, 152 x 193 cm	1
Order unit: PAK = Pack of 4 sets	




Contents	Pieces
Anesthesia drape, 183 x 300 cm	1
U-Drape, 193 x 304 cm	1
Cover drape, 102 x 150 cm	1
3-Pocket tray, 24 x 25 x 5 cm	1
Suture collector, 4 pockets	1
Preparation swab ball, middle, size 2, X-ray marker	5
Gauze swabs, plumsize, round, 20 x 20 cm, X-ray marker	10
Gauze compress, 16-fold (ply), 10 x 10 cm, X-ray marker	20
Disposable scalpel size 11	1
Askina® Elast fixation bandage, 4 m x 12 cm	1
Askina® Soft sterile, 7,5 x 5 cm	8
Askina® Soft sterile, 9 x 10 cm	1
Askina <sup>®</sup> universal bandage 5 m x 8 cm	1

Contents	Pieces
Askina <sup>®</sup> universal bandage 5 m x 10 cm	2
Suture Safil® violet 3/0 (2) 6 x 45 cm	1
Suture Safil® violet 0 (3.5) 70 cm, HR30	1
Suture Monosyn <sup>®</sup> undyed 4/0 (1.5) 45 cm, DS19	1
Suture Monosyn <sup>®</sup> undyed 4/0 (1.5) 70 cm, DS19	1
SUSI® needleholder Crile Wood, straight, 200 mm	1
SUSI® tissue forceps Adson, fine, 120 mm	1
SUSI® dissecting forceps Adson, fine, 120 mm	1
SUSI® Mosquito forceps, straight, 125 mm	3
SUSI® Mosquito forceps, curved, 125 mm	1
SUSI <sup>®</sup> dissecting scissors Metzenbaum, curved, 180 mm	1
SUSI® self-retaining retractor Weitlaner,	
3 x 4 teeth, 130 mm	1
SUSI <sup>®</sup> retractor Roux, 170 mm	1

Contents	Pieces
Abdominal drape, 6 ply, white, 45 x 45 cm, X-ray marker	5
Venostrip Vein exstirpations-set	1
Operation gown, Size L	3
Tablebag, 75 x 141 cm	1
Surgical gloves, size 8, powder free	1
Cover drape, 152 x 193 cm	1
Order unit: PAK = Pack of 4 sets	







#### SUSI® helps you in your daily work with many positive features and advantages

#### Safety

SUSI® reduces the risk of cross-contamination through surgical instruments.

#### Quality assurance

- With SUSI\* you always have a new instrument which ensures 100% functionality for the intended use
- SUSI<sup>®</sup> will be supplied sterile packed and are sterilised in a validated process

#### **Economical Aspects**

- SUSI® provides a clear cost allocation according to DRGs
- SUSI® leads to a reduction of processing costs in the hospital

#### Material

SUSI® are made from a high performance polymer. The special composition of this high-grade polymer ensures an outstanding quality and provides optimal mechanical properties for the intended function of each instrument. This material ensures a safe, easy and environmentally friendly removal through incineration in the regular hospital waste removal system.

#### Packaging

SUSI<sup>®</sup> are sterile packed in single peel packs. The single peel packs are presented in dispenser boxes containing an indicated quantity of units as described in this leaflet.





AN915SU	SUSI <sup>®</sup> Tubing Clamp
Dimension	135 mm, 5 <sup>5</sup> /16"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BF046SU	SUSI <sup>®</sup> Sponge Forceps
Name	GROSS
Description	straight, with universal ratchet
Dimension	200 mm, 8"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BF059SU	SUSI <sup>®</sup> Sponge Forceps
Description	curved, with universal ratchet
Dimension	265 mm, 10 <sup>1</sup> /2"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BF047SU	SUSI <sup>®</sup> Sponge Forceps
Name	GROSS
Description	curved, with universal ratchet
Dimension	200 mm, 8"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BC110SU	SUSI <sup>®</sup> Fine Scissors
Description	sh/sh, straight
Dimension	110 mm, 4 ¾"
Sales Unit	20 pieces, sterile, single packed in dispenser box




BC324SU	SUSI <sup>®</sup> Surgical Scissors
Description	sh/bl, straight
Dimension	145 mm, 5 ¾"
Sales Unit	20 pieces, sterile, single packed in dispenser box



1/1

BC606SU	SUSI <sup>®</sup> Dissecting Scissors
Name	METZENBAUM
Description	curved
Dimension	180 mm, 7"
Sales Unit	20 pieces, sterile, single packed in dispenser box

BC790SU	SUSI Umbilical Cord Scissors Blue	
BC791SU	SUSI Umbilical Cord Scissors Rose	
Description	curved	
Dimension	160 mm, 6 ¼ "	
Sales Unit	20 pieces, sterile, single packed in dispenser box	



BD222SU	SUSI® Dissecting Forceps
Name	ADSON
Description	w/o teeth, fine
Dimension	120 mm, 4 ¾
Sales Unit	20 pieces, sterile, single packed in dispenser box

BD512SU	SUSI <sup>®</sup> Tissue Forceps
Name	ADSON
Description	1 x 2 teeth, fine
Dimension	120 mm, 4 ¾"
Sales Unit	20 pieces, sterile, single packed in dispenser box






BD047SU	SUSI <sup>®</sup> Dissecting Forceps
Description	w/o teeth
Dimension	145 mm, 5 ¾"
Sales Unit	20 pieces, sterile, single packed in dispenser box

BD557SU	SUSI <sup>®</sup> Tissue Forceps
Description	1 x 2 teeth
Dimension	145 mm, 5 ¾"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BD239SU	SUSI® Dissecting Forceps	
Name	CUSHING	
Description	serrated	
Dimension	200 mm, 8"	
Sales Unit	20 pieces, sterile, single packed in dispenser box	

BD671SU	SUSI® Tissue Forceps
Name	WAUGHS
Description	1 x 2 teeth
Dimension	200 mm, 8"
Sales Unit	20 pieces, sterile, single packed in dispenser box



GK609SU	SUSI <sup>®</sup> Bipolar Coagulation Forceps	GK290SU -	GK292SU	SUSI <sup>®</sup> Bipolar Cable
Description	straight, w/o cable	Dimension	3.7 m	
Dimension	200 mm, 8" tip width 1.5 mm	Pin Distance	22 mm 28.6 mm	GK290SU (Dolly, EMC, ERBE-International) GK291SU (Aesculap, BOWA, EMC, Lamidey, Söring, Valleylab)
Sales Unit	20 pieces, sterile, single packed in dispenser box	Sales Unit	38 mm	GK292SU (Codman CMC II, CMC III) 10 pieces, sterile, single packed in dispenser box



BH110SU	SUSI <sup>®</sup> Mosquito Forceps
Name	HALSTED
Description	straight
Dimension	125 mm, 4 <sup>7</sup> / <sub>8</sub> "
Sales Unit	20 pieces, sterile, single packed in dispenser box



BH111SU	SUSI <sup>®</sup> Mosquito Forceps
Name	HALSTED
Description	curved
Dimension	125 mm, 4 <sup>7</sup> / <sub>8</sub> "
Sales Unit	20 pieces, sterile, single packed in dispenser box



BH196SU	SUSI <sup>®</sup> Haemostatic Forceps
Name	BIRKETT
Description	straight
Dimension	185 mm, 7 ¼"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BH197SU	SUSI <sup>®</sup> Haemostatic Forceps
Name	BIRKETT
Description	curved
Dimension	185 mm, 7 ¼"
Sales Unit	20 pieces, sterile, single packed in dispenser box




BH424SU	SUSI <sup>®</sup> Haemostatic Forceps
Name	PEAN
Description	straight
Dimension	145 mm, 5 ¾"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BM235SU	SUSI° Needle Holder
Name	MAYO-HEGAR
Dimension	150 mm, 6"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BM226SU	SUSI° Needle Holder
Name	CRILE WOOD
Dimension	200 mm, 8"
Sales Unit	20 pieces, sterile, single packed in dispenser box



BV067SU	SUSI® Self-Retaining Retractor
Name	WEITLANER
Description	semi-sharp, 3 x 4 teeth
Dimension	130 mm, 5 ¼"
Sales Unit	20 pieces, sterile, single packed in dispenser box

BT028SU	SUSI <sup>®</sup> Retractor
Name	ROUX
Description	blade sizes 36 x 20 mm & 20 x 15 mm
Dimension	170 mm, 6 ³/₄"
Sales Unit	20 pieces, sterile, single packed in dispenser box



FK850SU	SUSI <sup>®</sup> Curette
Description	sharp, round workpiece 7 mm oval workpiece 9 x 6 mm
Dimension	215 mm, 8 ½"
Sales Unit	20 pieces, sterile, single packed in dispenser box


<b>OM760SU</b>	SUSI® Tonsil Snare
Name	EVES
Description	with 2 wire loops
Dimension	285 mm, 11"
Sales Unit	10 pieces, sterile, single packed in dispenser box

PACK PER CASE	Product Code	Description	Page
	5061000	Suture Set 1	4
	5061010	Surgical Dressing Set 1	4
	5061020	Surgical Dressing Set 2	5
	5061030	Surgical Dressing Set 3	5
	5061040	Suture Removal 1	6
	5061050	Suture Removal	6
	5062000	Thorax Drainage Set, small	7
	5062010	Thorax Drainage Set, large	8
	5062020	Thorax Drainage Set	9
	5063000	SUSI <sup>®</sup> Port Set	10
	5063010	SUSI <sup>®</sup> Cardiac Pacemaker Set	11
	5063020	SUSI <sup>®</sup> Varicose Vein Set	12
CLINIC LINE	AN915SU	SUSI <sup>®</sup> Tubing Clamp	16
	BC110SU	SUSI <sup>®</sup> Fine Scissors	18
	BC324SU	SUSI <sup>®</sup> Surgical Scissors	18
	BC606SU	SUSI* Dissecting Scissors	19
	BD047SU	SUSI* Dissecting Forceps	21
	BD222SU	SUSI® Dissecting Forceps	20
	BD239SU	SUSI® Dissecting Forceps	22

### CLINIC LINE

Product Code	Description	Page
BD512SU	SUSI <sup>®</sup> Tissue Forceps	20
BD557SU	SUSI <sup>®</sup> Tissue Forceps	21
BD671SU	SUSI <sup>®</sup> Tissue Forceps	22
BF046SU	SUSI° Sponge Forceps	16
BF047SU	SUSI® Sponge Forceps	17
BF059SU	SUSI® Sponge Forceps	17
BH110SU	SUSI® Mosquito Forceps	24
BH111SU	SUSI® Mosquito Forceps	24
BH196SU	SUSI® Haemostatic Forceps	25
BH197SU	SUSI® Haemostatic Forceps	25
BH424SU	SUSI® Haemostatic Forceps	26
BM226SU	SUSI° Needle Holder	27
BM235SU	SUSI° Needle Holder	27
BT028SU	SUSI® Retractor	28
BV067SU	SUSI® Self-Retaining Retractor	28
FK850SU	SUSI <sup>®</sup> Curette	29
GK290SU GK291SU GK292SU	SUSI® Bipolar Cable	23
GK609SU	SUSI® Bipolar Coagulation Forceps	23
0M760SU	SUSI <sup>®</sup> Tonsil Snare	29





The main product trademark 'Aesculap' and the product trademark 'SUSI' are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. C82002

0314/0.6/7

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany Phone +49 (0) 74 61 95-0 | Fax +49 (0) 74 61 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company